


The following can be printed and used as reference and for recording foot measurements, which are then entered into LVP using the FOOTC data entry page. **Measurements should be entered in mm.** Measurements are most accurately taken using a digital caliper but any suitable method may be used. If you are interested in obtaining a suitable digital caliper contact us at support@limbvolumes.org for further information.

Ly	_____
Lxy	_____
W _x	_____
W _z	_____
W _y	_____
L12	_____
L8	_____
W12	_____
W8	_____
W4	_____
H _x	_____
H _z	_____
H _y	_____



For ease of data entry and consistency, enter data into the FOOTC page in millimeters (mm) only

**This page can be printed and used
as a Data Recording Sheet for
Your foot measurements**